

122 Monterey Avenue • Pacific Grove • California 93950 • (831) 375-0883 • (831) 375-0885 Fax Email: mdc@martincohen.com

CLEANING INSTRUCTIONS

THE MDC INTELLIGENT STIRRUPSTM ARE A PIECE OF MACHINERY WORKING IN A DIRTY ENVIRONMENT. ONCE A YEAR OR AS NEEDED, A NON GREASY, SPRAY LUBRICANT SHOULD BE APPLIED TO THE ROTATING TOP AND EXPOSED SCREW AND THEN WIPED COMPLETELY DRY. THE TWO OUTSIDE SET SCREWS CAN BE REMOVED AND THE INTERIOR OF THE STIRRUP CAN BE CLEANED AS WELL. THE SET SCREWS CAN THEN BE REPLACED, HOWEVER IF THEY ARE OVERTIGHTENED, THE MECHANISM WILL NOT TURN. ONE DROP OF BLUE, MEDIUM STRENGTH LOCTITE ADHESIVE WILL INSURE THEIR SECURE PLACEMENT. USERS SHOULD CHECK THEIR STIRRUPS REGULARLY AND KEEP THEM CLEAN AND FREE OF DEBRIS TO GUARANTEE PROPER WEAR AND MOVEMENT.

Thanks for your business and I hope you enjoy your new MDC Intelligent Stirrups,

MDC Corporation